


The Best Choice of **PIV** Engine in the World

# CACTUS 3.0

CACTUS is representative of multi-purpose **PIV** engine, especially oriented for the ultimate user convenience.



 **IIT** Co., Ltd.

# CACTUS 3.0

## IIT Co., Ltd.

Image Information Technology Co., Ltd.(IIT) was established in 1997. IIT is a research-based venture enterprise, specified in the field of image measurements and related software developments. Its final goal is to create technological power and vision for the full satisfaction of users. Main product is Particle Image Velocimetry(PIV) and related operating software engine PIV is the most advanced measurement technique ever developed to acquire multi-dimensional, simultaneous, multi-points, on-line, and real-time velocity information in all kind of flow fields which usually accompany complicated phenomena. PIV is an effective image measurement system to identify the accurate spacial coordinates of discrete tracer particles or particle neighbouring groups during two short time intervals. It always guarantees quantitative velocity vector information encountered in every flow fields of industrial and research purposes. In Korea, IIT has successfully developed the CACTUS series PIV software engines by support of the well-experienced and highly-motivated research team. In the future, IIT will concentrate its all innovative energy to pioneer engineering software development in the field of 4-D image measurement, and its comprehensive post-processing such as virtual reality animation based on Internet and multi-media environment.

## Outline of PIV

PIV (Particle Image Velocimetry) is the latest measurement technology which affords effective measurement of quantitative velocity vectors at multi-points and simultaneously from the visualized flow field by exploiting both the maximum CPU load of computer and modern image processing techniques. Adopting this technique, users can easily obtain the whole-field primitive velocity information of the time-dependent flow fields continuously and multi-dimensionally. Furthermore, the recorded flow image will be processed later at any convenient time without data loss, to satisfy reproducibility of the experimental data. PIV can be used as a complementary method with the point-wise measurement techniques such as LDV or hot-wire anemometry to secure local measurement accuracy. Nowadays, PIV is recognized as sole alternative to the Computational Fluid Dynamics(CFD) which furnishes various advantages impossible in the conventional experimental world.

## Applications of PIV

- Unsteady Turbulence Flow
- Aerodynamics (Vehicle, Aircraft, High-Speed Train etc.)
- Internal Flow of Engine
- Internal Flow of Heat Exchangers
- Internal Flow of Turbomachinery
- Wind Engineering
- Multiphase Flow
- Building Architecture and Environmental Flow
- Ocean Engineering
- Civil Engineering
- Bio-mechanics Flow
- Micro-Machine Flow
- Others



# Features of CACTUS 3.0

- Velocity Measuring Range:
  - 0,1mm/s ~ 400m/s(5W Argon-Ion)
  - Laser+Gated I.I. Camera, 100mmx100mm)
- Spatial Resolution:
  - 640x480 Pixel(Standard)
  - 2Kx2K(Max.)
- Sampling rate:
  - NTSC Video Signal : 60Hz
  - High-Speed Camera : Arbitrary
- Algorithm:
  - Non-FFT Cross-correlation PIV & Binary PTV
- Uncertainty < 0,32%
- High-End Performance
- Fast Processing Time:
  - 6.5 seconds for 40x40 grids on
  - Pentium-III 500MHz(interrogation
  - window size is 35x35 pixel and
  - maximum displacement is 8 pixel)
- Host Computer:
  - PC with built-in image grabber(DT series)
- O/S:
  - Win9X, Windows2000 or NT
- Large Processing Image Capacity:
  - max. RAM memory or extended
  - allowable hard disk memory
- Wide illumination sources:
  - Continuous Ar-Ion Laser
  - Nd-Yag Dual Pulse Laser
  - General White Lights
- Various Input Sources:
  - Analog CCD Camera
  - Digital Camera
  - Hi-8mm Camcorder
  - High Speed Video Camera
  - Gated Image Intensifier CCD
  - High Resolution Digital CCD(Up to 2K x 2K)
- 2-D Surface Flow Measurement at Arbitrary Camera Angle
- Multi-vision Measurement:
  - 3 CCD Camera (4:1 Aspect Ratio)
  - Mosaic Editing- Max. 15 Frames(H x V=5x3)
  - Support of Cinematic Animation

# CACTUS 3.0

## Merits of CACTUS 3.0

- Automatic on-line creation of AVI files for animation
  - Complete support for animation post-processing
  - Non-compressed base
- Velocity identification with high speed and real-time
- Complete arrangement of post-processing routines
  - High quality printout
  - Direct application of the original graphic software
- Minimum error occurrence
- High accuracy interpolation
- Support of user-friendly menu
- Programming tool : latest DELPHI visual language
- Continued version-up by module-type menu configuration

## Benefits of CACTUS 3.0 User


- Exclusive Membership of PIV research group
- Participation of PIV annual technical seminars arranged by IIT
- Support of optimized PIV configuration to meet user needs
- Provision of latest informations on PIV accessories and visualization equipments
- Support of 4-dimensional PIV and virtual reality-based post-processing in the future
- Technical support of complicated customer flows

## CACTUS 3.0 System



Items	Specification	Remark
PIV Software	- Animation support	- Fundamental specification of CACTUS 3.0 - Academic version available with major features
Image Grabber	- DT series	
Hi-8mm Camcorder	- SVHS	
Host Computer	- IBM compatible	- Compatible with DT board
High-speed Camera & Accessories	- 240 ~ 600 pps - B&W or RGB color	- Good choice of animation for high-speed flow
Argon-Ion Laser	- 4~6W(Water-cooling) - 300~500mW(Air-cooling)	- Continuous laser - Adequate for flow visualization
Nd-Yag Pulse Laser	- 30 ~ 300mJ/pulse	- Synchronization with pulse generator
Sheet-light Projector	- Optical fiber(10m) - Adjustable angle and thickness	- Generation of 2-Dimensional sheet-light(connected to laser)
CCD Camera	- B&W - 768x493 pixel	- For image input - Connected with Hi-8mm camcorder
Digital Camera	- 1024x1024 pixel	- For high-resolution input
Monitor	- 14 ~ 20 inch	- High resolution monitor
AOM	- 100MHz, 1 μs	- For continuous Ar-Ion Laser - Pulse Generator connected
Gated Image Intensifier Camera	- Gate time : 100ns - Max. Gain : $8,0 \times 10^4 \text{lm/m}^2/\text{lx}$	- Solution of the illumination problem for high-speed flow measurement, - Synchronization with pulse generator
Circulating water channel	- Velocity range: 0 ~ 0,5m/s, 300mmWx400mmH	- Supply of multi-purpose flow visualization
Particle	- Nylon12 - PVC - Expancel	- Tracer particle for water or air

# PIV System

## Arrangement of PIV System



## PIV Procedure







# Menu Description of CACTUS 3.0



CACTUS 3.0 Current Working Directory : C:/CACTUS/

File ImageProc PreProc ID PostProc View Window Tools Help



- File
- Load... Ctrl+O
- Close
- Close All
- Save As...
- Acquire...
- PassThru...
- Print... Ctrl+P
- Printer Setup...
- Preferences
  - General... Ctrl+Z
  - Change CACTUS Directory... Ctrl+D
- Exit





- ImageProc
- Zoom...
- Line Profile...
- Histogram...
- Filtering...
- Surface...



- PreProc
- Aspect Ratio Correction...
- Rotate & Scale Images...
- Transform OA Images...
- Make Background...
- Regress Log...
- Rotate Raw BGD...
- Coordinate Save...
- Obstacle Setting...



- ID
- Batch Mode...
- Identification...




# Menu Description of CACTUS 3.0


CACTUS 3.0 Current Working Directory : C:/CACTUS/

File ImageProc PreProc ID PostProc View Window Tools Help


- PostProc**
- MK Coord...
  - Auto Error Removing...
  - Manual Error Removing...
  - ReAllocation on Grid...
  - Convert To Real Unit...
  - Special Vector Average...
  - Averaging Vectors...
  - Compose Vectors...
  - u,v Component Extract...
  - Kinetic Energy...
  - Turbulence Intensity...
  - Vorticity...



- View**
- ID9 Vector...
  - Mkc Vector...
  - Aut Vector...
  - Man Vector...
  - Int Vector...
  - Vel Vector(Simple)...
  - Vel Vector...
  - Vel Vector(Lines)...
  - Contour (Lined & Filled)...
  - Surface...
  - Flow Lines...
  - Frequency Analysis...
  - X-Y Graph...



- Tools**
- Odd/Even Separation...
  - Sequential Rows...
  - Pathline...
  - Eval. Threshold...
  - Extract Raw8 from AVI...
  - Convert Bmp vs. Raw...
  - Ascii to Vel...
  - Binary to Ascii...
  - Wjf to Wmf...
  - File Attribute Change...
  - Rename ROK...
  - Palette Generation...
  - Create Various Pictures...
  - Create AVI clip...
  - Compose Bitmaps...
  - Labeling Images...





- Window**
- Tile Shift+Ctrl+F5
  - Cascade Shift+Ctrl+F4
  - Arrange Icons
  - ✓ 1 640x480-08 [OSA21.RAW]

- Help**
- Help Topics... F1
  - IT Home Page...
  - About...

# Application of CACTUS 3.0




## PIV Standard Project




## Refrigerator Evaporator Room


Original Image



Pathlines




Kinetic Energy




## Tundish


Velocity Vectors



Quasi-streamlines




Kinetic Energy




# Application of CACTUS 3.0

## Cylinder Wake Karman Vortex


Pathlines




Velocity Vectors



Streaklines




Quasi-streamlines




## Cleaning by Supersonic wave


Pathlines



Velocity Vectors




Quasi-streamlines




## Diffuser Pump Application of Gated I.I. Camera


Pump Geometry




Velocity Vectors



Kinetic Energy



Turbulence Intensity




# Application of CACTUS 3.0

## Measurement Conditions

- Algorithm : Cross-Correlation PIV
- Particle : CO<sub>2</sub>
- Laser : 300mW Ar-Ion Laser
- Input : I.I. Gated Camera
- Pulse Generator : 1  $\mu$ s Res.
- Pulse width : 6  $\mu$ s
- Pulse Interval : 6  $\mu$ s

Unit m/s	Min. Velocity	Max. Velocity	Mean Velocity
Instantaneous velocity field	36,5	149,6	82,1
Mean velocity field	49,0	98,6	75,2

## Arrangement of PIV System



## CO<sub>2</sub> Nozzle for Semi-conductor Cleaner


(A) Instantaneous

(B) Time Averaged


Original Image




Velocity Vectors




Velocity Profile




Path line




Velocity Vectors



Velocity Profile




Kinetic Energy




Kinetic Energy



Vorticity




Kinetic Energy



Kinetic Energy



Vorticity



## Surface Flow PIV

### (A) Groyne not used


Original Image




Transformed Image




Pathlines




Velocity Vectors



Velocity Profiles




Instantaneous Quasi-streamlines



### (B) Groyne used


Original Image




Transformed Image




Pathlines




Velocity Vectors



Velocity Profiles



Instantaneous Quasi-streamlines




# Application of CACTUS 3.0


## Pitching Aerofoil Multi-Vision PIV

Aerofoil : NACA0018, Attack Angle :  $18^\circ \sim 48^\circ$


Original Image




Instantaneous Velocity Vectors




Instantaneous Kinetic Energy




Instantaneous Vorticity




Instantaneous Turbulence Intensity




Instantaneous T.I. & Velocity




Instantaneous Turbulence K.E.




Instantaneous Reynolds Stress(Txx)




Instantaneous Reynolds Stress(Tyy)




Instantaneous Reynolds Stress(Txy)




Instantaneous Quasi-streamlines




Instantaneous Pressure Distribution




Velocity Profiles( $\Delta t=0.667$ )



Pathlines(grid :  $101 \times 61$ )



Streaklines(grid :  $41 \times 31$ )




# Animation of Pitching Aerofoil Multi-Vision PIV

Velocity Vectors

Quasi-streamlines

Pathlines



manufacturer

**IIT** Image Information Technology Co., Ltd.

Main office Dongjin INTEC Bldg, 1144-10, Choryang3-dong, Dong-ku, Pusan, 601-013, KOREA  
 TEL: 82-51-466-4601 FAX: 82-51-403-0381  
 LAB Flow Informatics Lab (N106), Korea Maritime University, Dongsam-dong1, Youngdo-ku, Pusan, 606-791, KOREA  
 TEL: 82-51-403-1213~4 FAX: 82-51-403-0381  
<http://www.iitpiv.com> e-mail: [info@iitpiv.com](mailto:info@iitpiv.com)



Chairperson, WWEC2009  
 Organizing Committee



Korea Maritime University

**Young-Ho LEE**

Professor / Ph.D.  
 Vice President, Korea Wind Energy Association (KWEA)

**Division of Mechanical Engineering**

1 Dongsam-dong, Youngdo-ku, Busan, 606-791, KOREA  
 TEL : (051)410-4293  
 FAX : (051)403-0381  
 Mobile : 010-3862-4293  
 E-mail : [lyh@hhu.ac.kr](mailto:lyh@hhu.ac.kr)  
<http://www.pivlab.net>